

# AUTOFLOW



CO<sub>2</sub> INCUBATORS



Best Products. Best Performance. Best Protection.

# The Perfect Choice for Your Unique Application



|                     | Model                     | Electronic Control System | Gas Control System  | Temperature Control System           | Humidity Control System                          | HEPA Filtration System | Configuration System                        |
|---------------------|---------------------------|---------------------------|---|--------------------------------------|--|------------------------|---|
| <b>IR AutoFlow</b>  | NU-8500, NU-8700          | Microprocessor            | CO <sub>2</sub> -Infrared   | Water-Jacketed                       | Convection <sup>6</sup>                          | Air/Gas <sup>4</sup>   | Single <sup>3</sup><br>Over-Under<br>Double |
| <b>US AutoFlow</b>  | NU-4750, NU-4850, NU-4950 | Microprocessor            | CO <sub>2</sub> -Infrared <sup>1</sup><br>O <sub>2</sub> Fuel Cell <sup>9</sup> | Water-Jacketed                       | Convection <sup>6</sup><br>Sensor <sup>2,8</sup> | Air/Gas <sup>4</sup>   | Single                                      |
| <b>DH AutoFlow</b>  | NU-5500                   | Microprocessor            | CO <sub>2</sub> -Infrared   | Direct Heat                          | Convection <sup>6</sup>                          | Air/Gas <sup>7</sup>   | Single <sup>3</sup>                         |
| <b>DHD AutoFlow</b> | NU-5510                   | Microprocessor            | CO <sub>2</sub> -Infrared   | Direct Heat with Sterilization Cycle | Convection <sup>6</sup>                          | Air/Gas <sup>7</sup>   | Single <sup>3</sup>                         |
| <b>TC PureCell</b>  | NU-5100                   | Microprocessor            | CO <sub>2</sub> -Thermo-Conductivity  | Direct Heat                          | Convection <sup>6</sup>                          | Air/Gas <sup>7</sup>   | Single <sup>3</sup>                         |


1] 2 or 3 gas control system  
 2] Monitors actual chamber humidity. Alarms in low humidity condition indicating that water reservoir requires filling  
 3] Right or left hand door swing – factory ordered

4] Internal recirculation system to ISO Class 5 clean air within 30 minutes  
 5] Right or left hand door swing – field reversible  
 6] Humidity pan (relative humidity to 95%)

7] Internal recirculation system to class 100 within 1 minute  
 8] Humidity control w / NU-4850, NU-4950  
 9] O<sub>2</sub> / N<sub>2</sub> control w / NU-4950  
**+ options**



Registered to: ISO 9001:2000  
 ISO 13485:2003



## Controlled Environments for Optimum Growth Conditions

- **The Sensitivity and Accuracy of Gas Control**

The AutoFlow incubators use a microprocessor-based, non-dispersive, digital solid-state infrared [IR] CO<sub>2</sub> sensor. This advanced design provides a very stable, drift-free output requiring less frequent calibrations.

- **The Precision of Electronic Microprocessor Controls**

AutoFlow incubators are run by a programmable microcontroller that samples input sensors, sets control outputs, and makes necessary corrections to the chamber environment. AutoFlow incubators provide quick and responsive recovery times.



- **The Long-Term Reliability and Protection of Automated and Redundant Back-Up Systems**

AutoFlow incubators provide dual sensor probes where appropriate, solid-state analyzers, multiple HEPA filtration systems, and the status alarms / outputs create a dependable, long-lasting scientific instrument.

- **The Flexibility of Available Options and Add-On Expansion Built into the System**

Additional shelves, RS-232 communication output, inner Lexan® doors, multi-signal analog outputs, tank switches, chart recorders, and moveable platforms are just a few of the available options and accessories.



### Constant Contamination Control

- Unique HEPA Filter System
  - Filters Fresh Air Into the Pump
  - Filters CO<sub>2</sub> Gas from Tanks or House Systems
  - Filters Internal Chamber Continuously to ISO Class 5 Air Quality
  - Filters IR CO<sub>2</sub> Sensor Sample Air Returning to Chamber
  - Filters Air Between Air Pump and Inner Chamber
- Interior Chamber is Constantly Maintained at Positive Pressure Similar to an ISO Class 5 Cleanroom
- Dual Sterilization Cycle (DHD AutoFlow NU-5510 Only)
  - 95°C Humidified Decontamination
  - 145°C Dry Sterilization Cycle





**NU-8700 / NU-8500**

# IR AutoFlow 8000 Series CO<sub>2</sub> Water-Jacketed Incubators

| Specifications               | Model          | Chamber Volume<br>(Ft. <sup>3</sup> / Liters) | Electrical*        | Chamber Dimensions<br>(W x H x D) | Exterior Dimensions<br>(W x H x D)                   | Weight                        |
|------------------------------|----------------|---|--------------------|-----------------------------------|--|-------------------------------|
|                              | <b>NU-8500</b> | 6.64 / 188                                    | 115 VAC / 60 Hz    | 21.5 x 25.5 x 21 in.              | 25.5 x 40.5 x 27 in.                                 | <b>Dry:</b> 218 lbs. / 99 kg  |
| <b>D:</b> 100 VAC / 50-60 Hz |                |   | 546 x 647 x 535 mm | 648 x 1029 x 686 mm               | <b>Full H<sub>2</sub>O Jacket:</b> 385 lbs. / 175 kg |                               |
| <b>E:</b> 230 VAC / 50-60 Hz |                |   |                    |                                   | <b>Shipping:</b> 287 lbs. / 130 kg                   |                               |
| <b>G:</b> 220 VAC / 50-60 Hz |                |   |                    |                                   |  |                               |
|                              | <b>NU-8700</b> | 13.28 / 376                                   | 115 VAC / 60 Hz    | 21.5 x 25.5 x 21 in.              | 25.5 x 73.5 x 27 in.                                 | <b>Dry:</b> 492 lbs. / 224 kg |
| <b>D:</b> 100 VAC / 50-60 Hz |                |   | 546 x 647 x 535 mm | 648 x 1867 x 686 mm               | <b>Full H<sub>2</sub>O Jacket:</b> 826 lbs. / 375 kg |                               |
| <b>E:</b> 230 VAC / 50-60 Hz |                |   |                    |                                   | <b>Shipping:</b> 561 lbs. / 255 kg                   |                               |
| <b>G:</b> 220 VAC / 50-60 Hz |                |   |                    |                                   |  |                               |

\*Specify NU-8500 / 8700, NU-8500D / 8700D, NU-8500E / 8700E, or NU-8500G / 8700G

## Features

### Standard Features

- 100% Stainless Steel Covered Interior Chamber
- 20 gal. / 75.7 L Capacity Water Jacket
- Removable Stainless Steel Shelves (4) [21 maximum] 19.25 in.<sup>2</sup> [489 mm<sup>2</sup>]
- Remote Alarm Output Contacts
- CO<sub>2</sub> Sample Port
- Water Level Sensor
- Adjustable Leg Levelers
- Access Port
- Fill Port and Drain Valve
- NuAire Incubator Control Electronics
- Full Size Water Pan (6.75 L)
- 8 ft. / 2.5 m Electrical Cord

### Optional Features

- Automatic CO<sub>2</sub> Tank Switch
- Communications Interface (RS-232)
- Internal Coil for Chilled Water
- Multi-Signal Analog Outputs: 0-5 VDC, 0-10 VDC, 4-20 mA
- 4 Inner Lexan® Doors
- Humidity Display (RH)

### Temperature Control System

- Default Set Point:** 37°C
- Chamber Temperature Range:** 18°C to 55°C (5°C Above Ambient to 30°C Ambient Max)
- Chamber Temperature Uniformity:** ± 0.2°C at 37°C
- Temperature Sensitivity:** ± 0.0125°C

### Electrical Requirements

- Model NU-8500:** 115 VAC / 60 Hz  
**Startup Power:** 625 Watts  
**Running Power:** 250 Watts  
**Heat Rejected:** 14 BTU/min. [Other electrical voltages available.]
- Model NU-8700:** 115 VAC / 60 Hz  
**Startup Power:** 1250 Watts  
**Running Power:** 500 Watts  
**Heat Rejected:** 28 BTU/min. [Other electrical voltages available.]

### CO<sub>2</sub> Control System

- Default Set Point:** 5%
- CO<sub>2</sub> Range:** 0% to 20%
- CO<sub>2</sub> Accuracy:** ± 0.1%
- CO<sub>2</sub> Recovery:** 5.0± 0.2% Within 3.5 Minutes
- CO<sub>2</sub> Sensor Type:** Non-dispersive Digital Infrared CO<sub>2</sub> Sensor

### Utility Connections

- Gas Connections:** .25" (6.3 mm) Tubing Connections
- Gas Input Pressure:** 20 PSIG (1.4 BAR) Input Pressures Maximum. Two-stage Gas Regulators Required.





# US AutoFlow 4000 Series CO<sub>2</sub> Water-Jacketed Incubators

| Specifications               | Model          | Chamber Volume<br>(Ft. <sup>3</sup> / Liters) | Electrical*        | Chamber Dimensions<br>(W x H x D)  | Exterior Dimensions<br>(W x H x D)                                  | Weight                        |
|------------------------------|----------------|---|--------------------|------------------------------------|---|-------------------------------|
|                              | <b>NU-4750</b> | 6.64 / 188                                    | 115 VAC / 60 Hz    | 21.5 x 25.5 x 21 in.               | 31.0 x 34.0 x 25.5 in.  | <b>Dry:</b> 222 lbs. / 101 kg |
| <b>D:</b> 100 VAC / 50-60 Hz |                |   | 546 x 647 x 535 mm | 787 x 864 x 648 mm                 | <b>Full H<sub>2</sub>O Jacket / RH system:</b><br>392 lbs. / 178 kg |                               |
| <b>E:</b> 230 VAC / 50-60 Hz |                |   |                    | <b>Shipping:</b> 291 lbs. / 132 kg |   |                               |
| <b>G:</b> 220 VAC / 50-60 Hz |                |   |                    |                                    |   |                               |
|                              | <b>NU-4850</b> | 6.64 / 188                                    | 115 VAC / 60 Hz    | 21.5 x 25.5 x 21 in.               | 31.0 x 34.0 x 25.5 in.  | <b>Dry:</b> 232 lbs. / 105 kg |
| <b>D:</b> 100 VAC / 50-60 Hz |                |   | 546 x 647 x 535 mm | 787 x 864 x 648 mm                 | <b>Full H<sub>2</sub>O Jacket / RH system:</b><br>412 lbs. / 188 kg |                               |
| <b>E:</b> 230 VAC / 50-60 Hz |                |   |                    | <b>Shipping:</b> 301 lbs. / 137 kg |   |                               |
| <b>G:</b> 220 VAC / 50-60 Hz |                |   |                    |                                    |   |                               |
|                              | <b>NU-4950</b> | 6.64 / 188                                    | 115 VAC / 60 Hz    | 21.5 x 25.5 x 21 in.               | 31.0 x 34.0 x 25.5 in.  | <b>Dry:</b> 232 lbs. / 105 kg |
| <b>D:</b> 100 VAC / 50-60 Hz |                |   | 546 x 647 x 535 mm | 787 x 864 x 648 mm                 | <b>Full H<sub>2</sub>O Jacket / RH system:</b><br>412 lbs. / 188 kg |                               |
| <b>E:</b> 230 VAC / 50-60 Hz |                |   |                    | <b>Shipping:</b> 301 lbs. / 137 kg |   |                               |
| <b>G:</b> 220 VAC / 50-60 Hz |                |   |                    |                                    |   |                               |

\*Specify NU-4750, 4850, 4950 / NU-4750D, 4850D, 4950D / NU-4750E, 4850E, 4950E / NU-4750G, 4850G, 4950G

## Features

### Standard Features

- 100% Stainless Steel Covered Interior Chamber
- 20 gal. / 75.7 L Capacity Water Jacket
- Removable Stainless Steel Shelves (4) [21 maximum] 19.25 in.<sup>2</sup> [489 mm<sup>2</sup>]
- Remote Alarm Output Contacts
- CO<sub>2</sub> Sample Port
- Water Level Sensor
- Adjustable Leg Levelers
- Access Port
- Fill Port and Drain Valve
- NuAire Incubator Control Electronics
- Full Size Water Pan (6.75 L)
- 8 ft. / 2.5 m Electrical Cord

### Optional Features

- Automatic CO<sub>2</sub> Tank Switch
- Communications Interface (RS-232)
- Internal Coil for Chilled Water
- Multi-Signal Analog Outputs: 0-5 VDC, 0-10 VDC, 4-20 mA
- 4 Inner Lexan® Doors
- Humidity Display (RH) and control (NU-4850 and NU-4950)
- O<sub>2</sub> Control System (NU-4950)

### Temperature Control System

- Default Set Point:** 37 °C
- Chamber Temperature Range:** 18 °C to 55 °C (5 °C Above Ambient to 30 °C Ambient Max)
- Chamber Temperature Uniformity:** ± 0.2 °C at 37 °C
- Temperature Sensitivity:** ± 0.0125 °C

### Electrical Requirements

- 115 VAC / 60 Hz
- Startup Power:** 625 watts
- Running Power:** 250 watts, 60 Hz
- Heat Rejected:** 14 BTU/min. [Other electrical voltages available.]

### Utility Connections

- Gas Connections:** .25" (6.3 mm) Tubing Connections
- Gas Input Pressure:** 20 PSIG (1.4 BAR) Input Pressures Maximum. Two-stage Gas Regulators Required.

### CO<sub>2</sub> Control Systems

- CO<sub>2</sub> (NU-4750 / 4850 / 4950)**
- Default Set Point:** 5%
- CO<sub>2</sub> Range:** 0 to 20%
- CO<sub>2</sub> Accuracy:** ± 0.1%
- CO<sub>2</sub> Recovery:** 5.0 ± 0.2% Within 3.5 Minutes
- CO<sub>2</sub> Sensor Type:** Non-dispersive Digital Infrared CO<sub>2</sub> Sensor

### RH (NU-4850 / 4950)

- Default Set-Point:** 90%
- RH Range:** 5% Above Ambient to 95%
- RH Accuracy:** ± 3%
- RH Recovery:** 95% ± 3%; Ten Minutes

### O<sub>2</sub> (NU-4950)

- Default Set-Point:** 21%
- O<sub>2</sub> Range:** 2 to 21%
- O<sub>2</sub> Accuracy:** ± 1.0%
- O<sub>2</sub> Recovery:** 5% ± 2%; Ten Minutes Average





**NU-5500**

# DH AutoFlow 5500 CO<sub>2</sub> Direct Heat Incubator

| Specifications | Model   | Chamber Volume<br>(Ft. <sup>3</sup> / Liters) | Electrical*           | Chamber Dimensions<br>(W x H x D) | Exterior Dimensions<br>(W x H x D) | Weight              |
|----------------|---------|---|-----------------------|-----------------------------------|------------------------------------|---------------------|
|                | NU-5500 | 6.64 / 188                                    | 115 VAC / 60 Hz       | 21.5 x 25.5 x 21.0 in.            | 25.5 x 39.5 x 26 in.               | 208 lbs. / 94.35 kg |
|                |         |   | E: 230 VAC / 50-60 Hz | 546 x 648 x 533 mm                | 648 x 1003 x 660 mm                |                     |
|                |         |   | G: 220 VAC / 50-60 Hz |                                   |                                    |                     |

\*Specify NU-5500 / NU-5500E / NU-5500G

## Features

### Standard Features

- 100% Stainless Steel Covered Foil Heaters Attached to Chamber Outer Walls
- Removable Stainless Steel Shelves (4) [17 maximum] 19.25 in.<sup>2</sup> [489 mm<sup>2</sup>]
- Remote Alarm Output Contacts
- CO<sub>2</sub> Sample Port
- Adjustable Leg Levelers
- Access Port
- Digital Infrared CO<sub>2</sub> Gas Sensor
- Water Pan (6.75 L)
- Microcomputer-Based Controls
- Replaceable HEPA Filter
- Replaceable Fan Wheel
- 8 ft. / 2.5 m Electrical Cord

### Optional Features

- Automatic CO<sub>2</sub> Tank Switch
- Communications Interface (RS-232)
- Multi-Signal Analog Outputs: 0-5 VDC, 0-10 VDC, 4-20 mA
- Additional Shelves
- Surge Protector

### Temperature Control System

- Temperature Range:** 18 °C to 55 °C (5 °C Above Ambient to 30 °C Ambient)
- Temperature Sensitivity:** ± .125 °C
- Temperature Uniformity:** ± 0.3 °C at 37 °C
- Temperature Accuracy:** ± 0.1 °C
- Temperature Recovery:** 0.3 °C / min. (Average)
- Temperature Display Resolution:** 0.1 °C

### Electrical Requirements

- 115 VAC / 60 Hz
- Startup Power:** 345 watts
- Running Power:** 175 watts
- Heat Rejected:** 10 BTU/min. [Other electrical voltages available.]

### Utility Connections

- Gas Connections:** .25" (6.3 mm) Tubing Connections
- Gas Input Pressure:** 20 PSIG (1.4 BAR) Input Pressures Maximum. Two-stage Gas Regulators Required.

### CO<sub>2</sub> Control Systems

- CO<sub>2</sub> Range:** 0.1 to 20%
- CO<sub>2</sub> Accuracy:** ± 0.1%
- CO<sub>2</sub> Recovery:** 5.0% ± 0.2% In Four Minutes
- CO<sub>2</sub> Display Resolution:** ± 0.1%
- CO<sub>2</sub> Sensor Type:** Non-dispersive Digital Infrared CO<sub>2</sub> Sensor



**[A]** A 99.98% efficient HEPA filter maintains Constant Contamination Control (C3) in the growth chamber.

**[B]** Dual temperature sensor probes, mounted just under the inflow to the HEPA filter controlling the chamber environment.





NU-5510

# DHD AutoFlow 5510 CO<sub>2</sub> Direct Heat Incubator with Decontamination Cycle

| Specifications | Model          | Chamber Volume<br>(Ft. <sup>3</sup> / Liters) | Electrical*   | Chamber Dimensions<br>(W x H x D)            | Exterior Dimensions<br>(W x H x D)          | Weight              |
|----------------|----------------|---|---|--|---|---------------------|
|                | <b>NU-5510</b> | 6.65 / 188.51                                 | 115 VAC / 60 Hz<br>E: 230 VAC / 50-60 Hz<br>G: 220 VAC / 50-60 Hz | 21.5 x 25.5 x 21.0 in.<br>546 x 648 x 533 mm | 25.5 x 39.5 x 26 in.<br>648 x 1003 x 660 mm | 208 lbs. / 94.35 kg |

\*Specify NU-5510 / NU-5510E / NU-5510G

## Features

### Standard Features

- 100% Stainless Steel Covered Interior Chamber
- Foil Heaters Attached to Chamber Outer Walls
- Removable Stainless Steel Shelves (4) [17 maximum] 19.25 in.<sup>2</sup> [489 mm<sup>2</sup>]
- Remote Alarm Output Contacts
- CO<sub>2</sub> Sample Port
- Adjustable Leg Levelers
- Access Port
- Digital Infrared CO<sub>2</sub> Gas Sensor
- Water Pan (6.75 L)
- Microcomputer-Based Controls
- Replaceable HEPA Filter
- Replaceable Fan Wheel
- 8 ft. / 2.5 m Electrical Cord
- Dual Decontamination Cycles**
  - 95°C Wet Decontamination Cycle
  - 145°C Dry Decontamination Cycle

### Optional Features

- Automatic CO<sub>2</sub> Tank Switch
- Communications Interface (RS-232)
- Multi-Signal Analog Outputs: 0-5 VDC, 0-10 VDC, 4-20 mA
- Additional Shelves
- Surge Protector

### Temperature Control System

- Temperature Range:** 18°C to 55°C (5°C Above Ambient to 30°C Ambient Max)
- Temperature Sensitivity:** ± .125°C
- Temperature Uniformity:** ± 0.3°C at 37°C
- Temperature Accuracy:** ± 0.1°C
- Temperature Recovery:** 0.3°C / min. (Average)
- Temperature Display Resolution:** 0.1°C

### Electrical Requirements

- 115 VAC, 50/60 Hz
- Startup Power: 345 watts
- Running Power: 175 watts
- Heat Rejected: 10 BTU/min.
- Decon. Cycle Power: 1150 watts [Other electrical voltages available.]

### Utility Connections

- Gas Connections: .25" (6.3 mm) Tubing Connections
- Gas Input Pressure: 20 PSIG (1.4 BAR) Input Pressures Maximum. Two-stage Gas Regulators Required.

### CO<sub>2</sub> Control System

- CO<sub>2</sub> Range:** 0.1 to 20%
- CO<sub>2</sub> Accuracy:** ± 0.1%
- CO<sub>2</sub> Recovery:** 5.0% ± 0.2% in Four Minutes.
- CO<sub>2</sub> Display Resolution:** ± 0.1%
- CO<sub>2</sub> Sensor Type:** Non-dispersive Digital Infrared CO<sub>2</sub> Sensor

### Decontamination Cycles

- 95°C Humidified Cycle**
- 95°C Humidified Decontamination Cycle Eradicates Contaminating Agents
- Advanced Chamber Design Allows for Shorter Heat Up and Cool Down Cycles; Complete Cycle Takes 14 Hours to Run
- Remote CO<sub>2</sub> IR Sensor Does Not Need to be Removed or Recalibrated
- 145°C Dry Cycle Features**
- 145°C High Temperature Dry Decontamination Cycle Eradicates Contaminating Agents
- Complete Cycle Takes 10 Hours to Run
- Remote CO<sub>2</sub> IR Sensor Does Not Need to be Removed or Recalibrated



[A]

[A] Covered interior corners prevent contaminants from becoming trapped in the chamber.





**NU-5100**

# TC PureCell 5100 CO<sub>2</sub> Direct Heat Incubator

| Specifications | Model          | Chamber Volume<br>(Ft. <sup>3</sup> / Liters) | Electrical*   | Chamber Dimensions<br>(W x H x D)            | Exterior Dimensions<br>(W x H x D)           | Weight                                 |
|----------------|----------------|---|---|--|--|--|
|                | <b>NU-5100</b> | 4.4 / 125                                     | 115 VAC / 60 Hz<br>E: 230 VAC / 50-60 Hz<br>G: 220 VAC / 50-60 Hz | 18.0 x 22.0 x 19.0 in.<br>457 x 559 x 533 mm | 23.0 x 30.0 x 23.4 in.<br>584 x 762 x 595 mm | 135 lbs. / 62.3 kg                     |
|                |                |   |   |  |  | *Specify NU-5100 / NU-5100E / NU-5100G |

## Features

### Standard Features

- 100% Stainless Steel Interior Chamber
- Foil Heaters Attached to Chamber Outer Walls
- Removable SST Shelves (3), Maximum Shelf Capacity (16)
- Remote Alarm Output Contacts
- CO<sub>2</sub> Sample Port
- Adjustable Leg Levelers
- Access Port
- Water Pan (4.8 liter Capacity)
- Microcomputer-Based Controls
- Replaceable HEPA Filter
- Replaceable Fan Wheel
- 8 ft. / 2.5 m Electrical Cord

### Optional Features

- Automatic CO<sub>2</sub> Tank Switch
- Communications Interface (RS-232)
- Multi-Signal Analog Outputs: 0-5 VDC, 0-10 VDC, 4-20 mA
- Additional Shelves
- Surge Protector

### Temperature Control System

- Temperature Range:** 18 °C to 55 °C (5 °C Above Ambient to 30 °C Ambient Max)
- Temperature Uniformity:** ± 0.5 °C at 37 °C
- Temperature Accuracy:** ± 0.25 °C
- Temperature Display Resolution:** 0.1 °C

### Electrical Requirements

- 115 VAC, 50/60 Hz
- Startup Power:** 130 watts
- Running Power:** 72 watts
- Heat Rejected:** 4.1 BTU/min. [Other electrical voltages available.]

### Utility Connections

- Gas Connections:** .25" (6.3 mm) Tubing Connections
- Gas Input Pressure:** 20 PSIG (1.4 BAR) Input Pressures Maximum. Two-stage Gas Regulators Required.

### CO<sub>2</sub> Control System

- CO<sub>2</sub> Range:** 0.1 to 20%
- CO<sub>2</sub> Accuracy:** ± 0.25%
- CO<sub>2</sub> Display Resolution:** ± 0.1%
- CO<sub>2</sub> Sensor Type:** Digital Thermo-conductivity



## Distributor Information

### To Locate a NuAire Distributor

Please visit NuAire's Web site at: [www.nuaire.com](http://www.nuaire.com) and select the 'Sales Team' menu item. Information can be found there for both United States and international distributors.



### Other NuAire Products



**LABGARD**

Biological Safety Cabinets



**GLACIER**

ULT Freezers



**AIREGARD**

Laminar Airflow Products



For more information contact NuAire, Inc. at 1.800.328.3352 or [www.nuaire.com](http://www.nuaire.com)



**NuAire, Inc.**

2100 Fernbrook Lane

Plymouth, MN 55447

U.S.A.

**1.800.328.3352**

**[www.nuaire.com](http://www.nuaire.com)**